



Panoramic Digital System

# Point 200 HD

Only Panoramic Model

Extra-Oral Imaging System





PointNix · HALLIM · OSLAM IMPLANT




# Point 200 HD

## Specification

### Point 200 HD

**Only Panoramic**

- The best efficient model among 'Only Panoramic device'
- High resolution Pano image by SMART HD Sensor

| Device      | Pano | CT Upgradable | Ceph Upgradable |
|-------------|------|---------------|-----------------|
| Point 200HD | 0    | X             | X               |


## Function

|                     |                                  |              |               |
|---------------------|----------------------------------|--------------|---------------|
| Panoramic Scan Time | 11sec with HD* mode              | Tube Current | 4 ~ 10mA      |
| Patient Position    | Standing / Wheel-Chair ccessible | Focal Spot   | 0.5mm X 0.5mm |
| Tube Voltage        | 50 ~ 90kvp                       | Data Bit     | 14 bit        |

\* HD mode: image enhancement function

## | Point 200 HD

### Adoption SMART HD Sensor




Easily show clear and high quality image.  
When the SMART HD Sensor absorbed X-ray, the each cell converts electronic signal without any spread or overlap.


#### Excellent rate of X-ray absorption

CCD - CMOS absorbs 35% of X-ray, but it needs more X-ray. However, SMART HD Sensor absorbs 100% of X-ray, and it shows high quality image with low X-ray dosage.

### Panoramic HD (High Definition) Image Quality with AMD (Automatic Management Dosage) solution



Panoramic image without HD mode




Panoramic image with HD mode

### CDX-View i version(Mobile Viewer)



Samsung GALAXY Note 10.1



iPad

# Dimension Information

## Point 200 HD



Model: Point 200 HD  
Weight: 150kg



PointNix · HALLIM · OSLAM IMPLANT

3F, #285, Digital-ro, Guro-gu, Seoul, Korea Tel : +82-2-854-8315(ext. 309) Fax : +82-2-839-4973  
E-mail : info@pointnix.com Web : www.pointnix.com, www.hallim.com, www.oslamimplant.com